

Soluscope Serie XS

Compact automated solution for small* endoscopes**



Soluscope

Rapid ISO 15883-1 & 4 validated cycle with full traceability

Standardise reprocessing to reduce the risk of cross-contamination

The automated and repeatable process of the *Soluscope Serie XS* is validated according to the strictest requirements, the ISO 15883-1 & 4 standards, to provide a **clean** and **disinfected** endoscope.

Ensure traceability of reprocessing activities

The Soluscope Serie XS includes a barcode reader, printer and local data storage and is compatible with the EP-Digital suite. It enables electronic traceability, real-time monitoring of your endoscope inventory, microbiological sampling and service records.

This allows you to comply with local guidelines by ensuring full traceability, thus contributing to reduce the risk associated with healthcare-related infections.







SOLUSCOPE PAA¹

- · Peracetic acid-based disinfectant
- Broad spectrum antimicrobial activity

SOLUSCOPE ECL Enzymatic Cleaner¹

- Enzymatic detergent
- Tested cleaning efficacy against biofilm and proteins²

Ergonomic solution designed to remove the need for laborious manual reprocessing



Streamline operator onboarding

With a large touchscreen and intuitive graphical interface, the *Soluscope Serie XS* is easy-to-use, helping staff working with efficiency.



Contribute to improve staff satisfaction

The specific height of the machine enables a comfortable and ergonomic position for users when operating.



Safety-focused design to prevent any direct contact with chemicals

The Lock & key system and integrated canulas on containers enable to prevent any direct contact with chemicals. Additionally, the Soluscope Serie XS operates without requiring a fume extractor or specialized room ventilation, thereby preventing staff from exposure to chemical fumes



Simple to install and adapted to small rooms

Compact, mobile and with minimal installation requirements, the *Soluscope Serie XS* can be easily installed in small facilities.



Automated solution to help give time back to staff

Enhance efficiency and productivity, with minimal staff intervention

In just 19 minutes, it provides a cleaned and disinfected endoscope for the next examination, so as to support the pace of the daily operating program, even with minimal staff availability.

Lower the total cost of ownership by protecting endoscopes

The endoscope is fully immersed with a rotating flow so that all the surfaces of the endoscope are properly reprocessed.

The continuous leak testing controls the integrity of the shaft during the cycle, protecting the endoscopes in case of damage of the outer surface.



Contribute to overall sustainability

Ecolab places sustainable innovation at the heart of our medical solutions to support our customers in their commitments to the overall sustainability. Every Soluscope equipment and chemical product are designed to minimize its environmental impact while delivering optimal performance.



REDUCED PACKAGING

Concentrated chemistries reduce consumption by 50%³ while maintaining high-level performance.



ENERGY SAVINGS


As a result of a responsible industrial approach, the energy consumption of the *Soluscope Serie XS* has been reduced by 33%³.



RESPONSIBLE MANUFACTURING

All of our equipment and chemistry are responsibly designed and manufactured in Europe, focusing on quality and reliability.





HELPS TO DELIVER













Ecolab delivers a complete solution to improve the reprocessing of the flexible endoscopes at every step of the process. To help protect patients' safety and reduce the risk of contamination, each product is designed for simplicity, ergonomics and efficiency.

Thanks to this combined offering, each endoscope reprocessing department benefits from improved user wellbeing, operational efficiency and sustainability, for enhanced patient care. Lean on Ecolab's expert knowledge, service and training and enhance your experience with a complete solution.

FIND SIMPLICITY IN SAFETY AND MASTERY IN EFFICIENCY

- *Flexible ENT endoscopes and bronchoscopes
- **Semi-critical, heat-sensitive flexible endoscopes with or without channels



SOLUSCOPE SAS

(an Ecolab Company) 100, Rue du Fauge Z.I Les Paluds - 13400 Aubagne - France

www.soluscope.com

Distributed by:

ECOLAB EUROPE GMBH

Hofwiesenstrasse 349 CH-8050 Zürich, Switzerland

www.ecolab.com







D03194A